

UMI-Gen: A UMI-based read simulator for variant calling evaluation in paired-end sequencing NGS libraries

Vincent Sater, Pierre-Julien Viailly, Thierry Lecroq, Philippe Ruminy, Caroline Bérard, Élise Prieur-Gaston, Fabrice Jardin

▶ To cite this version:

Vincent Sater, Pierre-Julien Viailly, Thierry Lecroq, Philippe Ruminy, Caroline Bérard, et al.. UMI-Gen: A UMI-based read simulator for variant calling evaluation in paired-end sequencing NGS libraries. Computational and Structural Biotechnology Journal, 2020, 18, pp.2270-2280. 10.1016/j.csbj.2020.08.011. hal-03210438

HAL Id: hal-03210438 https://cnrs.hal.science/hal-03210438

Submitted on 5 Sep 2022

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



Distributed under a Creative Commons Attribution - NonCommercial 4.0 International License

UMI-Gen: a UMI-based read simulator for variant calling evaluation in paired-end sequencing NGS libraries

Vincent Sater^{a,c,d,e}, Pierre-Julien Viailly^{b,c,d}, Thierry Lecroq^a, Philippe Ruminy^{b,c}, Caroline Bérard^a, Élise Prieur-Gaston^a, Fabrice Jardin^{b,c}

^a University of Rouen Normandy UNIROUEN, LITIS EA 4108, 76000 Rouen, France
 ^b Department of Pathology, Centre Henri Becquerel, 76000 Rouen, France
 ^c INSERM U1245, University of Rouen Normandy UNIROUEN, 76000 Rouen, France
 ^d These authors have contributed equally to this work
 ^e To whom correspondence should be addressed

Abstract

Motivation: With Next Generation Sequencing becoming more affordable every year, NGS technologies asserted themselves as the fastest and most reliable way to detect Single Nucleotide Variants (SNV) and Copy Number Variations (CNV) in cancer patients. These technologies can be used to sequence DNA at very high depths thus allowing to detect abnormalities in tumor cells with very low frequencies. Multiple variant callers are publicly available and are usually efficient at calling out variants. However, when frequencies begin to drop under 1%, the specificity of these tools suffers greatly as true variants at very low frequencies can be easily confused with sequencing or PCR artifacts. The recent use of Unique Molecular Identifiers (UMI) in NGS experiments has offered a way to accurately separate true variants from artifacts. UMI-based variant callers are slowly replacing raw-read based variant callers as the standard method for an accurate detection of variants at very low frequencies. However, benchmarking done in the tools publication are usually realized on real biological data in which real variants are not known, making it difficult to assess their accuracy.

Results: We present UMI-Gen, a UMI-based read simulator for targeted sequencing paired-end data. UMI-Gen generates reference reads covering the targeted regions at a user customizable depth. After that, using a number of control files, it estimates the background error rate at each position and then modifies the generated reads to mimic real biological data. Finally, it

Preprint submitted to Journal Name

July 30, 2020

will insert real variants in the reads from a list provided by the user. Availability: The entire pipeline is available at https://gitlab.com/vincentsater/umigen under MIT license. Contact: vincent.sater@gmail.com

Keywords: Sequence Analysis, UMI, Simulator, Variant Calling, NGS

1 1. Introduction

Nowadays, next generation sequencers such as Thermo Fisher or Illumina 2 have become the standard go-to method for DNA sequencing. Prior to se-3 quencing, DNA must be extracted and amplified by PCR in order to generate 4 enough fragments to cover the wanted amplicons. After amplification, the 5 sequencer handles the obtained fragments and generates their sequences in 6 form of reads. In most applications, especially ones that handle variant detection, the obtained reads must then be aligned to a reference genome in 8 order to be used effectively. Today, cancer diagnosis is a very active area of 9 research and one of its most important applications is the detection of Single 10 Nucleotide Variants (SNV) in tumor cells. In fact, each cancer type has a 11 specific profile of genetic mutations in specific genes. Therefore, establishing 12 a precise profile of variants in a cancer patient allows to better understand the 13 cancer evolution and customize the treatment according to the established 14 profile. 15

Detecting and calling out variants in the aligned reads is done through a 16 variant calling analysis. Generally, variant calling tools can detect mutational 17 events such as substitutions, insertions and deletions very efficiently. How-18 ever, at very low variant allele frequencies (VAFs) (under 1%), it becomes 19 very challenging for raw-read-based variant callers to accurately call variants. 20 In fact, PCR amplification and the sequencing step can introduce errors in 21 the final reads. These errors are called artifacts and occur at very low VAFs 22 which can lead to the confusion between them and true low-frequency vari-23 ants. Multiple studies [1, 2, 3, 4, 5] have shown the effectiveness of using 24 Unique Molecular Identifiers as a way to filter out PCR and sequencing arti-25 facts. UMIs are short arbitrary oligonucleotide sequences that are attached 26 to DNA fragments by ligation before the PCR amplification. By definition, 27 the UMI tags must be random sequences so each fragment can have a unique 28 short oligonucleotide sequence attached to it, giving each fragment a unique 29

sequence tag. During the amplification, the UMI tags are amplified with 30 their respective fragments. After sequencing, each UMI tag can be figured 31 out from the reads. The idea behind using UMI tags in NGS experiments 32 to filter out artifacts is explained in Figure 1. In fact, if a variant is a true 33 mutation, it means that it must have been present on the initial DNA frag-34 ment so when we tag the DNA fragment with a UMI, we are also tagging the 35 mutation. The fragments that result from the amplification of that mutated 36 DNA fragment must all be tagged by the same UMI tag and carry the same 37 mutation (Figure 1A). On the other hand, if the variant is a sequencing error, 38 it means that the initial DNA fragment did not have the mutation in the first 39 place and that it appeared later in the sequencing step. Therefore, during 40 the amplification step, all the fragments resulting from the amplification of 41 that DNA fragment should theoretically be tagged with the same UMI and 42 should not present the mutation. The mutation will be produced later on, 43 in the sequencing step, affecting only some reads and not all of them, thus 44 creating discrepancies in the same UMI group (Figure 1B). 45

With the growing number of variant calling tools, it has become hard to 46 choose the right tool adapted to a certain experiment. Data simulation can 47 play an important role for testing different tools on a dataset that we have 48 control on, a control that we do not have on real biological data. At the 49 moment, many short read simulators exist such as IntSIM [6] that can simu-50 late somatic variants using HMM models trained on real sequencing genomes 51 and SVSR [7] that is specifically designed to simulate datasets with structural 52 variations and is compatible with multiple sequencing platforms. These tools 53 are publicly available for researchers and allow them to test their algorithms 54 on a simulated dataset in which variants are inserted at different frequencies 55 and at different positions. The usage of the read simulators enable having a 56 very accurate benchmarking of each variant calling tool ability. Surprisingly, 57 no simulation software exists at the moment that let users generate reads 58 with UMI tags. In this article, we present UMI-Gen, a UMI-based read sim-59 ulator that can be used not only to test raw-read based variant callers but 60 most importantly, UMI-based ones. UMI-Gen uses multiple real biological 61 samples to estimate background error rate and base quality scores at each 62 position. Then, it will introduce real variants in the final reads. To test our 63 tool, we used 6 control samples and show exactly how our algorithm estimate 64 the background error rate at each position. Then we give it a list of 15 vari-65 ants at different positions and at different frequencies to introduce them in 66 the final reads. Finally, we used 2 raw-read-based variant callers: SiNVICT 67



Figure 1: The difference between a true variant and an artifact from a UMI perspective. (A) A true variant is present on the DNA fragment so when the UMI tag 1 is added, it tags the fragment and the mutation as well. After amplification, all the fragments tagged with the UMI tag 1 carry the same mutation. (B) An artifact is not present on the DNA fragment but rather appears at the steps that follow the UMI introduction. That is why not all fragments with the same UMI tag 2 carry the same artifact.

[8] and OutLyzer [9] and two UMI-based variant callers: DeepSNVMiner
[10] and UMI-VarCal [11] in order to compare the 4 tools performance and
demonstrate that UMI-Gen correctly inserts the given variants at their respective positions and at the correct frequencies in a dataset that mimics
perfectly what is seen in biological samples.

73 2. Materials and methods

74 2.1. Software input

UMI-Gen requires a minimum of three parameters at execution: a list of control BAM/SAM samples, the BED file with the coordinates of the targeted genomic regions and a reference genome FASTA file with BWA index files. In fact UMI-Gen is designed to work on targeted sequencing data only thus a BED file is always required. UMI-Gen can also accept a fourth optional file under the PILEUP format. In fact, when running UMI-Gen on control samples, a PILEUP file is automatically produced. This file contains the A, C, G and T average counts at each position for all the control
samples. This file can be given to UMI-Gen at execution time and will allow
the software to reload the pileup generated during the last analysis instead
of regenerating it. This will allow the user to gain some significant time since
the pileup generation is the most time-consuming step.

87 2.1.1. Control samples

Control samples are BAM/SAM files that are obtained by sequencing healthy individuals and normally should not contain any somatic variant. UMI-Gen can accept input files in BAM and SAM formats. A pileup is performed on each sample and a final average pileup is generated from the counts of all control samples.

93 2.1.2. Variant file

This file contains a list of the variants the user wishes to insert in the simulated reads. These are the only variants that should be reported in the variant callers VCF file during variant calling benchmarks. The variant file is a Comma Separated Values (CSV) file that contains 2 columns: the first column contains the variant ID with the HGVS nomenclature and the second column being the variant's desired frequency. UMI-Gen will then go to each position and insert these variants in order to produce final reads.

101 2.2. Generating the final pileup

102 2.2.1. Pileup

The first step of the workflow (Figure 2) consists of generating the final 103 pileup. For each control sample, our pileup algorithm will count the occur-104 rences of each A, C, G and T. The counts will be stored for each position 105 of the BED file as well as the average quality of the position and its depth. 106 This is basically the same algorithm that is used by UMI-based variant caller 107 UMI-VarCal that is been reintegrated in this tool for its high efficiency in 108 treating reads with UMI tags. When all the pileups for all the control sam-109 ples are ready, they will be merged in a final pileup that contains the average 110 statistics (counts, depth and quality score) at each position based on the 111 observations on all control samples (Figure 2A). When the average pileup is 112 complete and ready, it will be automatically dumped as a PILEUP file that 113 contains all the calculated information on the set of control samples. If the 114 user wishes to generate simulated data based on the same BED file and the 115 same set of control samples, the dumped pileup can be used directly which 116



Figure 2: Background error estimation workflow. (A) The first step runs over every position in all control samples and counts the total occurrences of every A, C, G and T. It also stores the average base quality score for each position. (B) The second step's goal is to remove any suspected variant from the pileup as our objective is to estimate background error noise only. (C) In this step, the counts are converted to probabilities by dividing them by the depth for each position. (D) The final step consists of converting the base quality score of each position to the corresponding ASCII + 33 character.

allows the program to skip the pileup generation step and go directly to thevariant calling step, saving the user much significant time.

119 2.2.2. Variant calling

Even though the control samples are theoretically variant-free, SNP and undetected mutations could still be present in the files. These potential variants must be removed so they would not be present in the final reads. To do so, we used the same variant calling method implemented in UMI-VarCal to call out potential variants and remove them from the pileup. This step will produce what we call a filtered pileup (Figure 2B).

¹²⁶ 2.2.3. Background noise estimation

The background noise estimation step consists of calculating the frequency of observing an A/C/G/T at each position. Without the background errors, at each position the reference base should have a frequency of 1 while

the remaining three bases should be at 0. The total of the four frequen-130 cies must be equal to 1. However, we know that artifacts exist in our control 131 samples and these artifacts represents the background noise that we normally 132 encounter in a normal NGS experiment. Since our aim is to simulate reads 133 that are highly similar to those produced with real sequencing experiments, 134 UMI-Gen calculates the real base frequencies from the control samples at 135 each position. The frequencies will then be used as a probability matrix 136 when producing the final reads. When this step is complete, a probability 137 pileup is generated (Figure 2C). Insertions and deletions are not considered 138 during the background noise estimation and thus, are not present in the final 139 pileup as their occurrence has a much lower rate (~ 1000 times lower) than 140 that of substitutions) especially in second and third generation sequencers 141 [12]. Therefore, we judge that their inclusion is not worth complicating the 142 algorithm for. 143

144 2.2.4. Quality scores estimation

Our tool was developed on sequencing files produced by an Illumina se-145 quencer. In the FASTQ files produced by Illumina sequencers, quality scores 146 are encoded into a compact form, which uses only 1 byte per quality value 147 [13, 14]. The full table of encoding is available in Table S1. UMI-Gen is there-148 fore only compatible with sequencers that use the same encoding. UMI-Gen 140 calculates the average quality score for each position based on the qualities in 150 all control samples and then converts the quality score to the corresponding 151 ASCII character to be inserted in the final FASTQ file. This is the final step 152 of the pileup generation workflow and will produce the final pileup (Figure 153 2D). Moreover, UMI-Gen also models the base quality scores per position in 154 read on the control samples and introduces the estimation in the final reads. 155 Based on all the reads in the control samples, our tool will calculate a median 156 base quality score for each position in the reads to produce a quality per po-157 sition matrix. This matrix is then used at the end to recalibrate the quality 158 scores according to each base's position in the read. For example, this allows 159 UMI-Gen to mimic the loss of quality at the end of the reads when present. 160

161 2.3. Producing the reads

The main objective of UMI-Gen is to generate paired-end reads that mimic reads obtained from real life experiments. To do so, it starts exactly the way a real-life sequencing experiment starts: getting the DNA fragments. At the beginning, our tool will generate a number of initial sequences that



Figure 3: The difference between adding a true variant and adding an artifact in generated reads. (A) Adding an artifact is relatively easy as all the tool has to do is to modify the base at the wanted position without touching the read's UMI tag. (B) On the other hand, in order to add a true variant, the software must change the base at the wanted position on a set of reads. Then it will create a new UMI tag (UMI tag 3) and change the UMI tag of all the affected reads to UMI tag 3.

only present the reference base at each position. The user can explicitly 166 specify the desired length for all the reads at execution. It should be noted 167 that the algorithm will only create reads that will exactly align on the speci-168 fied positions from the BED file so off-target amplification is not considered. 169 Then, a UMI tag is attached to each initial sequence. Depending on the 170 amplification factor and the desired depth chosen by the user, the algorithm 171 will keep amplifying the initial sequences until the desired depth is reached at 172 all positions. In fact, at this step, default values for the amplification factor 173 and initial DNA fragments are automatically calculated in order to ensure 174 optimal performance of the tool. We do so by analyzing the depth chosen 175 by the user and the VAFs of the variants that he wishes to introduce. Using 176 these numbers, we calculate the minimum number of initial DNA fragments 177 needed for the true variant insertion. Even though this will ensure optimal 178 performance, the user is free to change these parameters as long as they 179 are mathematically allowed. Once we have the reference reads, the second 180

step consists of adding the background noise (refer to section 2.2.3) to these 181 reads (Figure 3A). Using the probability matrix calculated before, UMI-Gen 182 modifies the reads at each position for them to match the calculated prob-183 abilities. These modifications are done without changing the reads' UMI 184 so they mimic PCR and sequencing artifacts: they are false positives and 185 should not be called by variant callers. Finally, UMI-Gen parses the variant 186 file provided by the user in order to insert true mutations in the final reads. 187 The algorithm will go to each position, change the probability of the variant 188 to the corresponding frequency from the variant file. In this step, since UMI-180 Gen is adding a true variant, the UMI tags of the modified reads are also 190 modified in order to produce concordant UMI tags (Figure 3B). A concor-191 dant UMI tag is a UMI whose all reads carry the exact same mutation. Also, 192 since UMI-Gen generates paired-end data, when adding a mutation on one 193 read, the variant is automatically added to its mate (since we only generate 194 paired reads that always overlap). 195

196 2.4. Software output

Once all variants are inserted, UMI-Gen will generate the two FASTQ 197 files (R1 and R2). It will then call BWA [15] to do the alignment, a step that 198 will produce a BAM file. SAMtools [16] is finally called to create the BAM's 199 index file and convert the BAM into SAM. All five files are generated in the 200 desired output directory. In addition, UMI-Gen generates a binary PILEUP 201 file that corresponds to the dumped average pileup. This file can be used to 202 skip the pileup regeneration and load the pileup directly if the analysis was 203 already done on the same control samples. 204

205 2.5. Implementation

Launching UMI-Gen's workflow (Figure 4) is handled by a main Python 206 script that controls many Python3 modules. In order to achieve better over-207 all performance, Cython was used to compile all Python modules. UMI-Gen 208 requires for the tools BWA and SAMtools to be installed on the PC/server: 209 BWA is called for the alignment step and SAMtools for converting, sorting 210 and indexing the generated BAM files. Our tool can be executed through 211 a UNIX/Linux command line interface. In total, UMI-Gen can accept 20 212 parameters at execution. Managing these parameters allows the user to have 213 full control over his simulated data. A list of all the parameters and thresh-214 olds is available in Table S2. 215



Figure 4: UMI-Gen's workflow: Control samples are used to create a background noise frequency matrix and the user provides a CSV file with a list of the wanted variants. Using the FASTA and the BED files, UMI-Gen creates a first set of UMI-tagged reference reads. Artifacts are then inserted to mimic the sequencer's background noise. Finally, the tool uses the list provided by the user to insert variants at their exact locations.

216 3. Results

217 3.1. Control samples

A targeted sequencing panel was designed at the Centre Henri Becquerel 218 in Rouen (France) to search for specific mutations in the DNA of patients 219 suffering from Diffuse Large B cell Lymphoma (DLBCL). This panel of 76,630 220 bases is designed to identify genomic abnormalities within a list of 36 genes 221 that are most commonly impacted in this type of lymphoma. The panel 222 was specifically designed for QIAseq chemistry allowing UMI introduction 223 in the DNA fragments during the construction of the library. A list of the 224 genes used in the panel and their corresponding number of targeted regions 225 is provided in the supplementary Table S3. In order to test our tool's ability 226 to mimic and reproduce average sequencer background noise in the produced 227 sample, we randomly selected 6 samples from a very large number of patients 228 whose DNA were sequenced at the Centre Henri Becquerel. All six samples 229 are liquid biopsies with circulating cell-free DNA that was checked to be 230 adequate for sequencing. We preferred the use of liquid biopsies as these 231

samples usually contain a high number of very low frequency variants and
artifacts. Using such samples as control samples will produce simulated data
with a relatively high number of artifacts. This will allow us to have an
accurate estimate of the specificity of each tested variant caller.

Sample	Α	С	G	Т
Control 1	0	11	10	874
Control 2	0	1	7	843
Control 3	0	2	2	860
Control 4	0	6	9	965
Control 5	1	2	4	867
Control 6	3	2	2	880

Table 1: The A,C,G and T breakdown at position 2,493,165 of chromosome 1 for the six control samples.

Table 1 shows the exact counts of A, C, G, T for position 2,493,165 236 on chromosome 1 for each control sample. The fist control sample counts 237 (0,11,10,874), the second sample has (0,1,7,843), the third one has (0,2,2,860), 238 the fourth sample shows (1,6,9,965), the fifth one has (1,2,4,867) and the final 239 one counts (3,2,2,880). As explained in section 2.2.3, UMI-Gen will calculate 240 an average count for each base and then estimate its probability. In our 241 case and for this position, the obtained average count has 4 A, 24 C, 34 G 242 and 5289 T with a total count of 5351 bases. To obtain the probabilities 243 for this position on this chromosome, we simply divide each base count by 244 the total count of the 4 bases, obtaining the final probability vector (0.0007, 245 0.0045, 0.0064, 0.9884). If, for example, we wanted to produce a BAM file 246 with a depth of 3000x, this position would have 2 A, 14 C, 19 G and 2965 T. 247 The probability matrix mentioned in section 2.2.3 is basically the probability 248 vectors of each position of the panel, merged together. In our test and in 249 order to demonstrate our results, we simulated two artificial samples in which 250 we added the calculated background error noise. The first sample or Sample 251 1 has an average depth of 1000x (+/-15%) at each position) and Sample 252 2 has an average depth of 10,000x. To make sure that the artifacts were 253 correctly added to the reads, we used IGV (version 2.4.16) [17] to visualize 254 the reads. Figure 5 shows how the background error noise is properly and very 255 accurately added at position 2,493,165 of chromosome 1 with the probabilities 256 calculated from the 6 control samples above. 257



Figure 5: The A, C, G and T breakdown at the position 2,493,165 of the chromosome 1 in the produced samples: Sample 1 with the depth of 1000x (A) and Sample 2 with the depth of 10,000x (B).

258 3.2. Simulated data validation

In order to validate our simulated dataset, we compared it to the control 259 samples used to generate it. First, we compared the base quality scores 260 distribution in the reads. Figure 6A shows the variation of the median base 261 quality scores with the position of base in the read for the control samples. 262 We can clearly see that the median score is very high and very stable at the 263 start and all along the read's length (>34). However, a first drop in quality 264 is noted at position 138 and a second more considerable one at position 145. 265 In our simulated data, we chose an average length for the reads of about 110 266 bp so the longest read had a length of 127. We can see, in Figure 6B, how the 267 algorithm perfectly recreates the stability of the scores all along the simulated 268 reads. However, since the simulated reads did not have lengths > 135 bp, 269 we do not see that little drop at the end of the simulated reads. In fact, 270 to be sure that our quality score estimation works correctly, we simulated a 271 drop in quality at the position 85 and wanted to see if it will be inserted in 272 the simulated reads. Figure 7 shows how the simulated drop in quality (38) 273 \rightarrow 34) at position 85 was perfectly reproduced in the simulated data (36 \rightarrow 274 33). Another parameter we wanted to verify is the %GC variation between 275



Figure 6: The variation of the median base quality score with position in read in real samples (A) and in the simulated data (B).

the control and the simulated data. Figure 8 clearly shows how the median %GC of reads in the control data (Figure 8A - 56% GC) is nearly identical to that of the simulated reads (Figure 8B - 57% GC).

279 3.3. Inserted variants

Two different lists of mutations were created to go along with each simu-280 lated sample. The first list contains 11 substitution variants with frequencies 281 that go from 0.9 (90%) to 0.01 (1%), one deletion at 1% and one insertion 282 at 1%. This list is used to produce the simulated Sample 1 with a depth 283 of 1000x. The second list contains 13 substitution variants with frequencies 284 that go from 0.9 (90%) to 0.001 (0.1%), one deletion at 1% and one insertion 285 at 1%. This list is used to produce the simulated Sample 2 with a depth of 286 10,000x. Two very low frequency variants (frequency < 1%) were added to 287 the second list to test the variant insertion accuracy of UMI-Gen. In fact, 288



Figure 7: The variation of the median base quality score with position in read in real samples (A) and in the simulated data (B). A simulated drop in quality was simulated in scenario A and its reproduction in the simulated dataset (B).

very low frequency variants are the hardest to detect and should be system-289 atically used to rigorously test any variant caller. In order to verify that the 290 wanted variants were added at the exact locations with the correct frequen-291 cies, we used IGV to visualize the reads. Figure 9 shows the variants added 292 in both samples and Table 2 details the exact variants that we inserted at 293 the specific locations. Next generation sequencers have difficulties with accu-294 rately detecting variants in long homopolymer regions. Some variant callers 295 automatically filer out variants that occur in such regions and others do not. 296 In order to avoid any bias, we chose each variant's location carefully to make 297 sure that it is not inserted in a homopolymer region. Figure 10 demonstrates 298 that our tool is capable of accurately adding variants in the final reads at 299 the specified locations for both samples. 300



Figure 8: The repartition of the GC percentage in reads in the real data (A) and in the simulated data (B).

301 3.4. Variant detection

We tested the ability of four different variant callers to correctly detect the 302 true variants added in section 3.3 and filter out sequencing errors/artifacts 303 added in 3.1. We used SiNVICT and OutLyzer, two raw-read-based variant 304 callers specifically developed to detect low frequency variants and two UMI-305 based variant callers (DeepSNVMiner and UMI-VarCal) with a very low fre-306 quency detection threshold and that analyze UMI tags in order to produce 307 more accurate results. The four variant callers were tested on the two artifi-308 cial samples: Sample 1 that contains 13 known variants and a depth of 1000x 309 and Sample 2 that contains 15 known variants and a depth of 10,000x. Both 310 samples have a total of 76,630 sequenced positions which corresponds to the 311 size of the sequencing panel. 312

Table 3 and 4 detail the results of each tool for Sample 1 and 2 respectively. The total number of positives corresponds to the number of variants

Position	Reference allele	Variant allele	Frequency	Sample
2,488,101	G	А	0.9	S1 & S2
$2,\!489,\!200$	С	А	0.8	S1 & S2
2,491,260	А	G	0.7	S1 & S2
2,493,201	Т	А	0.6	S1 & S2
$2,\!494,\!300$	G	А	0.5	S1 & S2
23,885,600	С	А	0.4	S1 & S2
23,885,800	А	Т	0.3	S1 & S2
27,022,900	С	А	0.2	S1 & S2
27,023,200	С	А	0.1	S1 & S2
27,093,001	G	А	0.05	S1 & S2
27,100,350	С	А	0.01	S1 & S2
$27,\!106,\!500$	G	А	0.005	S2 only
117,057,400	Т	А	0.001	S2 only
120,458,000	С	CTA	0.1	S1 & S2
120,466,600	TGTC	Т	0.1	S1 & S2

Table 2: Detailed list of the inserted mutations. In this test, all mutations are inserted on chromosome 1.

found in the result VCF file. The total number of negatives is then calculated 315 by subtracting total positives from the total number of positions (76,630). 316 The four variant callers had comparable results between the two samples. 317 Starting with SiNVICT, it detected 241 variants in Sample 1 and 463 in 318 Sample 2 but with the same number of true positives. This corresponds to a 319 sensitivity of 61.5%/53.4% which is relatively acceptable and a specificity of 320 99.7%/99.4% on Sample 1/2. Moving on to OutLyzer, the tool detected 109 321 variants in Sample 1 and three times more variants in Sample 2 (342). Unfor-322 tunately, this corresponded to one more true positive, the rest being only false 323 positives. Outlyzer scored good sensitivities (> 80%) and excellent specifici-324 ties (99.9%/99.6%) on both samples. Concerning DeepSNVMiner, the tool 325 managed to detect all the inserted variants except the deletion in both sam-326 ples. The tool scored very high scores on sensitivity (92.3%/93.4%) as well 327 as specificity (99.95%/99.99%) for both datasets. Finally, UMI-VarCal was 328 able to achieve a perfect score (100%) in terms of sensitivity and specificity 329 on both samples detecting all the 13/15 variants in Sample 1/2 with no false 330 positives for both configurations. 331



Figure 9: Along with the reference genome FASTA file and the BED file, two different lists were used, one with 13 variants and the other with 15 variants to respectively produce the artificial samples Sample 1 and Sample 2.

332 3.5. Performance

In order to evaluate UMI-Gen's performance, we simulated four samples 333 with increasing depths: 500, 1000, 5000 and 10,000. For each simulated 334 sample, execution time and memory consumption were reported. The four 335 samples were simulated using the same six control samples. The first time 336 we run UMI-Gen, the pileup generation step is mandatory. The pileup gen-337 eration step only depends on the control samples and takes about 1.5 minute 338 per sample. The quality estimation step following the pileup is also essential 339 and takes on average 0.5 minute per sample. However, these 2 steps gener-340 at files that can be given directly to the program at the execution. This 341 means that for the other times the user wants to simulate data using the 342 same control samples, the pileup file and the quality matrix file can be used 343 directly allowing to save considerable time. Table 5 details the execution 344 time numbers and the memory needed to generate each sample. Generating 345 the FASTQ files takes only 1.57 minute for the 500x sample and uses only 346 1 GB of RAM. On the other side, 16.58 minutes are needed for a sample of 347 10,000x and memory consumption goes up to 5.1 GB. All these tests were 348 performed on a computer running Linux (Ubuntu 16.04) using only one core 349



Figure 10: The inserted mutations were correctly added to the reads with their exact locations at their corresponding frequencies. Here, we see four mutations: chr1:2491260A>G at 70%, chr1:27022900C>A at 20%, chr1:120458000C>CTA at 10% and chr1:27093001G>A at 5%.

CPU running at 2.20 GHz and equipped with 16 GB of RAM. All measurements were done three times and the average was used for the comparison.
After the FASTQ generation, BWA and SAMtools are called from within the tool to generate the corresponding BAM and SAM files.

354 4. Discussion

Tagging DNA fragments with UMI tags have proved itself as a very reli-355 able method to significantly reduce - if not completely remove - the number 356 of false positives upon variant calling. A huge number of variant callers 357 are publicly available at the moment but unfortunately, only 4 of them are 358 specifically developed to treat UMI tags in reads. For raw-read-based vari-359 ant callers, a lot of artificial read simulator exist and can satisfy everyone's 360 needs. However, to our knowledge, no tool is publicly available to simulate 361 artificial reads with UMI tags. Such tool is very important as it allows de-362 velopers to accurately test the specificity and the sensitivity of their variant 363 callers on artificial reads in which real variants are known instead of testing 364

Variant Caller	TP	FP	FN	Sensitivity (%)	Specificity (%)
Sinvict	8	233	5	61.5	99.7
OutLyzer	11	98	2	84.6	99.9
DeepSNVMiner	12	37	1	92.3	99.95
UMI-VarCal	13	0	0	100	100

Table 3: Variant calling results on Sample 1. Four variant callers were tested: SiNVICT, OutLyzer, DeepSNVMiner and UMI-VarCal and for each tool, True Positives (TP), False Positives (FP), False Negatives (FN), sensitivity and specificity are reported.

Variant Caller	TP	FP	FN	Sensitivity (%)	Specificity (%)
Sinvict	8	455	7	53.4	99.4
OutLyzer	12	330	3	80	99.6
DeepSNVMiner	14	2	1	93.4	99.99
UMI-VarCal	15	0	0	100	100

Table 4: Variant calling results on Sample 2. Four variant callers were tested: SiNVICT, OutLyzer, DeepSNVMiner and UMI-VarCal and for each tool, True Positives (TP), False Positives (FP), False Negatives (FN), sensitivity and specificity are reported.

them on biological samples whose mutational profile is completely or partiallyunknown.

Our main objective was to develop a UMI-based read simulator that is 367 fast, accurate and reliable. UMI-Gen is able to estimate the background 368 error noise of a given control dataset and then reproduce it accurately in the 369 produced reads. Doing so, it allows to mimic the sequencer's background 370 noise of a real sequencing experiment. We also showed that our simulator is 371 able to accurately insert variants if provided with a list of variants with exact 372 locations and their corresponding frequencies and produce reads that mimic 373 ones produced in real life experiments. In our tests, we were able to insert 374 mutations as low as 0.1% but theoretically, we can go as low as we want 375 provided that the depth of the produced sample is accordingly increased. 376

Moreover, in our variant caller comparison, SiNVICT did a decent job detecting the 8 of the added variants and went as low as 5%. Impressively, we judge the performance of OutLyzer as excellent as it detected 12 of the 15 variants (Sample 2) and showed a detection threshold of 0.5% which is very respectable. However, SiNVICT and OutLyzer being raw-read-based

Sample	Data Simulation (min)	FASTQ to BAM (s)	Ram Usage (GB)
500x	1.57	8	1.0
1000x	1.87	14	1.1
5000x	6.97	52	2.6
10,000x	16.58	99	5.1

Table 5: Performance analysis of UMI-Gen: the variation of execution time and memory consumption with the simulated data's depth.

variant callers, UMI tags were not treated in the reads and therefore, both 382 tools produced a high percentage of false positives. On the other hand, Deep-383 SNVM iner results were near perfect as expected from a decent UMI-based 384 variant caller detecting all variants except one in both scenarios with only a 385 couple of false positives. Finally, UMI-VarCal was successfully able to treat 386 UMI tags allowing it to filter out all false positives and only call out the 387 13 added variants in Sample 1 and all of the 15 in Sample 2. These results 388 demonstrate how the UMI-based variant calling approach is much more effi-389 cient and accurate than raw-read-based ones allowing to detect variant with 390 VAFs as low as 0.1% without sacrificing specificity. It also highlights the 391 need to the development and usage of UMI-based read simulators in order 392 to test these new algorithms. 393

³⁹⁴ 5. Conclusion

Here, we present UMI-Gen: a standalone UMI-based read simulator for 395 variant calling evaluation in paired-end sequencing NGS libraries. UMI-Gen 396 produces sequencing files (FASTQ, BAM and SAM) for an artificial sample 397 to be used for UMI-based variant calling testing purposes. By using a set of 398 control DNA samples, our tool is capable to accurately mimic the background 399 error noise of the sequencer and add it into the reads. After that, it can insert 400 specific mutations at specific locations and at very precise frequencies that 401 can go as low as 0.1% (and even lower). In our tests, all added artifacts were 402 correctly inserted in the reads, causing a high number of false positives in the 403 raw-read-based variant callers results. Also, all inserted true variants were 404 visualized with a genome visualizing tool (IGV) and were detected by at least 405 one of the four variant calling tools we tested. Finally, we note that UMI-406 Gen's filters and parameters (such as read length and UMI tag length) are 407

customizable which gives the user total control over his produced samples.
This level of customization allows the tool to be adequate for a high number
of research applications.

411 6. Funding

This work was partly funded by the Université de Rouen Normandie and Vincent Sater is funded by a PhD fellowship from the Région Normandie.

414 7. References

- M. W. Schmitt, S. R. Kennedy, J. J. Salk, E. J. Fox, J. B. Hiatt, L. A. Loeb, Detection of ultra-rare mutations by next-generation sequencing, Proc Natl Acad Sci U S A 109 (2012) 14508-14513. URL: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3437896/. doi:10.1073/pnas.1208715109.
- 418 [2] Y. Kukita, R. Matoba, J. Uchida, T. Hamakawa, Y. Doki, F. Imamura, K. Kato, High-fidelity target
 419 sequencing of individual molecules identified using barcode sequences: de novo detection and absolute
 420 quantitation of mutations in plasma cell-free DNA from cancer patients, DNA Res 22 (2015) 269-277.
 421 URL: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4535617/. doi:10.1093/dnares/dsv010.
- [3] A. M. Newman, A. F. Lovejoy, D. M. Klass, D. M. Kurtz, J. J. Chabon, F. Scherer, H. Stehr,
 C. L. Liu, S. V. Bratman, C. Say, L. Zhou, J. N. Carter, R. B. West, G. W. Sledge, J. B. Shrager,
 B. W. Loo, J. W. Neal, H. A. Wakelee, M. Diehn, A. A. Alizadeh, Integrated digital error suppression for improved detection of circulating tumor DNA, Nat Biotechnol 34 (2016) 547-555. URL:
 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4907374/. doi:10.1038/nbt.3520.
- 427 [4] A. L. Young, G. A. Challen, B. M. Birmann, T. E. Druley, Clonal haematopoiesis harbour 428 ing AML-associated mutations is ubiquitous in healthy adults, Nat Commun 7 (2016). URL:
 429 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4996934/. doi:10.1038/ncomms12484.
- [5] D. Z. Bar, M. F. Arlt, J. F. Brazier, W. E. Norris, S. E. Campbell, P. Chines, D. Larrieu, S. P.
 Jackson, F. S. Collins, T. W. Glover, L. B. Gordon, A novel somatic mutation achieves partial rescue in a child with Hutchinson-Gilford progeria syndrome, J Med Genet 54 (2017) 212-216. URL:
 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5384422/. doi:10.1136/jmedgenet-2016-104295.
- K. Yuan, J. Zhang, L. Yang, IntSIM: An Integrated Simulator of Next-Generation Sequencing Data, IEEE
 Transactions on Biomedical Engineering 64 (2017) 441–451. doi:10.1109/TBME.2016.2560939.
- [7] X. Yuan, M. Gao, J. Bai, J. Duan, SVSR: A Program to Simulate Structural Variations and Generate Sequencing Reads for Multiple Platforms, IEEE/ACM Transactions on Computational Biology and Bioinformatics 17 (2020) 1082-1091. doi:10.1109/TCBB.2018.2876527.
- (8) C. Kockan, F. Hach, I. Sarrafi, R. H. Bell, B. McConeghy, K. Beja, A. Haegert, A. W. Wy-att, S. V. Volik, K. N. Chi, C. C. Collins, S. C. Sahinalp, SiNVICT: ultra-sensitive detection of single nucleotide variants and indels in circulating tumour DNA, Bioinformatics 33 (2017) 26–34. doi:10.1093/bioinformatics/btw536.
- [9] E. Muller, N. Goardon, B. Brault, A. Rousselin, G. Paimparay, A. Legros, R. Fouillet,
 O. Bruet, A. Tranchant, F. Domin, C. San, C. Quesnelle, T. Frebourg, A. Ricou, S. Krieger,
 D. Vaur, L. Castera, OutLyzer: software for extracting low-allele-frequency tumor mutations
 from sequencing background noise in clinical practice, Oncotarget 7 (2016) 79485-79493. URL:
 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC5346729/. doi:10.18632/oncotarget.13103.

- [10] T. D. Andrews, Y. Jeelall, D. Talaulikar, C. C. Goodnow, M. A. Field, DeepSNVMiner: a sequence analysis tool to detect emergent, rare mutations in subsets of cell populations, PeerJ 4 (2016). URL: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4888318/. doi:10.7717/peerj.2074.
- [11] V. Sater, P.-J. Viailly, T. Lecroq, E. Prieur-Gaston, E. Bohers, M. Viennot, P. Ruminy, H. Dauchel,
 P. Vera, F. Jardin, UMI-VarCal: a new UMI-based variant caller that efficiently improves low-frequency
 variant detection in paired-end sequencing NGS libraries, Bioinformatics (Oxford, England) (2020).
 doi:10.1093/bioinformatics/btaa053.
- 455 [12] M. Schirmer, R. DAmore, U. Z. Ijaz, N. Hall, C. Quince, Illumina error pro456 files: resolving fine-scale variation in metagenomic sequencing data 17 (????). URL:
 457 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4787001/. doi:10.1186/s12859-016-0976-y.
- 458 [13] B. Ewing, P. Green, Base-Calling of Automated Sequencer Traces Using Phred. II. Error
 459 Probabilities, Genome Research 8 (1998) 186-194. URL: http://genome.cshlp.org/content/8/3/186.
 460 doi:10.1101/gr.8.3.186.
- [14] B. Ewing, L. Hillier, M. C. Wendl, P. Green, Base-Calling of Automated Sequencer
 Traces UsingPhred. I. Accuracy Assessment, Genome Research 8 (1998) 175-185. URL:
 http://genome.cshlp.org/content/8/3/175. doi:10.1101/gr.8.3.175.
- 464 [15] H. Li, R. Durbin, Fast and accurate short read alignment with burrows-wheeler transform,
 465 Bioinformatics 25 (2009) 1754-1760. URL: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2705234/.
 466 doi:10.1093/bioinformatics/btp324.
- 467 [16] H. Li, B. Handsaker, A. Wysoker, T. Fennell, J. Ruan, N. Homer, G. Marth, G. Abecasis, R. Durbin,
 468 The Sequence Alignment/Map format and SAMtools, Bioinformatics 25 (2009) 2078-2079. URL:
 469 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2723002/. doi:10.1093/bioinformatics/btp352.
- 470 [17] J. T. Robinson, H. Thorvaldsdóttir, W. Winckler, M. Guttman, E. S. Lander, G. Getz,
 471 J. P. Mesirov, Integrative Genomics Viewer, Nature biotechnology 29 (2011) 24-26. URL:
 472 https://www.ncbi.nlm.nih.gov/pmc/articles/PMC3346182/. doi:10.1038/nbt.1754.