

Complete sequences of two Paenibacillus sp. strains and one Lysinibacillus strain isolated from the environment on STAA medium highlight biotechnological potential

Sabrina A Attéré, Laurie C Piché, Laurent Intertaglia, Raphaël Lami, Steve J Charette, Antony T Vincent

▶ To cite this version:

Sabrina A Attéré, Laurie C Piché, Laurent Intertaglia, Raphaël Lami, Steve J Charette, et al.. Complete sequences of two Paenibacillus sp. strains and one Lysinibacillus strain isolated from the environment on STAA medium highlight biotechnological potential. Microbiology Resource Announcements, 2024, 13 (6), pp.00199-24. 10.1128/mra.00199-24. hal-04682272

HAL Id: hal-04682272 https://cnrs.hal.science/hal-04682272v1

Submitted on 30 Aug 2024

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.







3 | Environmental Microbiology | Announcement

Complete sequences of two *Paenibacillus* sp. strains and one *Lysinibacillus* strain isolated from the environment on STAA medium highlight biotechnological potential

Sabrina A. Attéré, 1,2,3 Laurie C. Piché, 1,4 Laurent Intertaglia, 5 Raphaël Lami, 6 Steve J. Charette, 1,2,3 Antony T. Vincent 1,4

AUTHOR AFFILIATIONS See affiliation list on p. 2.

ABSTRACT Streptomycin thallous acetate actidione medium is typically used to isolate *Brochothrix thermosphacta* bacteria from food. Using this medium, three bacterial strains were isolated from the environment. Genomic sequences demonstrated that these bacteria are of the genera *Lysinibacillus* and *Paenibacillus* and are of biotechnological interest.

KEYWORDS biotechnology, environmental microbiology, *Lysinibacillus*, *Paenibacillus*, secondary metabolites, STAA medium

Streptomycin thallous acetate actidione (STAA) agar is used to isolate *Brochothrix* thermosphacta bacteria (1), a meat spoilage agent (2). This medium has also been used to isolate soil bacteria, leading to the identification of *Brochothrix campestris* (3), the only other species in the *Brochothrix* genus apart from thermosphacta.

A sampling campaign was carried out in autumn of 2023 in Quebec City (Province of Quebec, Canada). Two 7.5-mL soil samples from different sites were each combined with 7.5 mL of water and left at room temperature for 6 h. An aliquot of each sample (1 mL) was spread on STAA agar medium and incubated at 25°C for 48 h. A water sample from a watering hole located in a wooded area was also plated on STAA agar medium and underwent incubation at 25°C for 48 h. White colonies were replated onto HIB agar and incubated at 25°C for 48 h. One isolate per sample (i.e., Y5S-7 [soil], Y5S-8 [soil], and Y5S-9 [water]), was frozen at -80°C in HIB liquid medium with 20% glycerol as a cryoprotectant.

Strains Y5S-7, Y5S-8, and Y5S-9 were thawed on HIB agar medium and incubated for 24-48 h at 25°C. The total DNA of the strains was extracted with QIAamp PowerFecal Pro DNA Kit (QIAGEN, Toronto, ON, Canada) following the manufacturer's instructions. DNA samples were sequenced by Plasmidsaurus (Eugene, OR, USA) using an Illumina NextSeq2000 (2 × 150 bp) and a Nanopore PromethION. The Illumina libraries were prepared using the SeqWell ExpressPlex 96 library prep kit, while the Nanopore libraries were prepared using v14 library prep chemistry without fragmentation or size selection, and sequenced on an R10.4.1 flow cell. Basecalling was carried out for Nanopore reads using Dorado version 7.1.4 (4) on super-accurate mode. Illumina reads were filtered using Fastp version 0.23.4 (5) with default parameters, while Nanopore reads were filtered using Filtlong version 0.2.1 (6) and retained the best 90% of reads above 1,000 bp or until only 500 Mbp remained. After filtration, 13.28 million, 5.84 million, and 9.03 million Illumina reads remained for Y5S-7, Y5S-8, and Y5S-9, respectively. A total of 86,759 (N50 = 6,329 bp), 78,181 (N50 = 6,804 bp), and 98,146 (N50 = 6,131 bp) Nanopore reads were retained for strains Y5S-7, Y5S-8, and Y5S-9, respectively. Read quality was verified with FastQC version 0.12.1 (7). A hybrid assembly was then performed for each data set using Unicycler version 0.5.0 (8). The sequences were linearized to the dnaA gene for chromosomes and repA for plasmids using Unicycler version 0.5.0 (8), and were annotated by the

Editor Irene L. G. Newton, Indiana University, Bloomington, Bloomington, Indiana, USA

Address correspondence to Antony T. Vincent, antony.vincent@fsaa.ulaval.ca.

The authors declare no conflict of interest.

See the funding table on p. 2.

Received 28 February 2024 Accepted 9 April 2024 Published 29 April 2024

Copyright © 2024 Attéré et al. This is an open-access article distributed under the terms of the Creative Commons Attribution 4.0 International license.

June 2024 Volume 13 Issue 6

TABLE 1 Genomic characteristics of Y5S-7, Y5S-8, and Y5S-9 strains

Genus	Species	Strain	Replicon	Length (bp)	CDSs	GC (%)	Coverage (×)	GenBank
Lysinibacillus	sp.	Y5S-8	Chromosome	4,672,069	4,442	37.78	194	CP145889
			pY5S8-1	4,614	6	41.46	46	CP145890
			pY5S8-2	2,217	2	41.5	80	CP145891
Paenibacillus	sp.	Y5S-9	Chromosome	6,924,530	5,927	45.94	195	CP145887
			pY5S9-1	4,690	5	44.43	63	CP145888
Paenibacillus	amylolyticus	Y5S-7	Chromosome	6,917,693	5,933	45.97	268	CP145892
			pY5S7-1	374,980	272	38.55	354	CP145893
			pY5S7-2	68,778	83	39.56	1,079	CP145894

NCBI using PGAP (9) on submission to GenBank. The assembly characteristics are shown in Table 1.

Taxonomic identification using the Type Strain Genome Server (10) suggests that Y5S-8 and Y5S-9 may be novel species within the *Lysinibacillus* and *Paenibacillus* genera, respectively, while Y5S-7 was identified as a *Paenibacillus amylolyticus* species. Annotation with antiSMASH (11) demonstrated that strains Y5S-7 and Y5S-9 have all the genes to produce bacillopaline, a metallophore (12), and polymyxin, a last-resort antibiotic (13) with biotechnological potential.

ACKNOWLEDGMENTS

This work was supported by the Natural Sciences and Engineering Research Council of Canada (Grant no. RGPIN-2022–03321) and the Samuel-de-Champlain program, joint between the FRQNT and the Consulat général de France of Québec.

AUTHOR AFFILIATIONS

¹Institut de Biologie Intégrative et des Systèmes, Pavillon Charles-Eugène-Marchand, Université Laval, Quebec City, Quebec, Canada

²Centre de Recherche de l'Institut Universitaire de Cardiologie et de Pneumologie de Québec, Hôpital Laval, Quebec City, Quebec, Canada

³Département de biochimie, de microbiologie et de bio-informatique, Faculté des Sciences et de Génie, Université Laval, Quebec City, Quebec, Canada

⁴Département des sciences animales, Faculté des Sciences de l'Agriculture et de l'Alimentation, rue de l'Agriculture, Université Laval, Quebec City, Quebec, Canada

⁵Sorbonne Université, CNRS, Bio2mar, Observatoire Océanologique de Banyuls-sur-Mer, Avenue Pierre Fabre, Banyuls-sur-Mer, France

⁶Sorbonne Université, CNRS, Laboratoire de Biodiversité et Biotechnologies Microbiennes, Observatoire Océanologique de Banyuls-sur-Mer, Avenue Pierre Fabre, Banyuls-sur-Mer, France

AUTHOR ORCIDs

Steve J. Charette https://orcid.org/0000-0002-0199-2852 Antony T. Vincent http://orcid.org/0000-0002-1338-5785

FUNDING

Funder	Grant(s)	Author(s)
Canadian Government Natural Sciences and Engineering Research Council of Canada (NSERC)	RGPIN-2022-03321	Antony T. Vincent
FRQ Fonds de recherche du Québec – Nature et		Steve J. Charette
technologies (FRQNT)		Antony T. Vincent
Consulat de France		Raphaël Lami

Downloaded from https://journals.asm.org/journal/mra on 30 August 2024 by 194.167.19.129.

DATA AVAILABILITY

The chromosomal and plasmid sequences of Y5S-7, Y5S-8, and Y5S-9 strains have been deposited in GenBank (Table 1). Illumina and Nanopore sequencing reads were deposited in the Sequence Read Archive (SRA) database: Y5S-7 (SRX23698269 and SRX23698272), Y5S-8 (SRX23698270 and SRX23698273), and Y5S-9 (SRX23698271 and SRX23698274).

REFERENCES

- Gardner GA. 1985. Streptomycin-thallous acetate-actidione (STAA) agar: a medium for the selective enumeration of *Brochothrix thermosphacta*. Int J Food Microbiol 2:69–70. https://doi.org/10.1016/0168-1605(85)90058-3
- Illikoud N, Rossero A, Chauvet R, Courcoux P, Pilet M-F, Charrier T, Jaffrès E, Zagorec M. 2019. Genotypic and phenotypic characterization of the food spoilage bacterium *Brochothrix thermosphacta*. Food Microbiol 81:22–31. https://doi.org/10.1016/j.fm.2018.01.015
- Talon R, Grimont PAD, Grimont F, Gasser F, Boeufgras JM. 1988. Brochothrix campestris sp. nov. Int J Syst Bacteriol38:99–102. https://doi. org/10.1099/00207713-38-1-99
- Oxford nanopore. 2023. Available from: https://github.com/nanopore-tech/dorado
- Chen S, Zhou Y, Chen Y, Gu J. 2018. fastp: an ultra-fast all-in-one FASTQ preprocessor. Bioinformatics 34:i884–i890. https://doi.org/10.1093/ bioinformatics/bty560
- Andrews S. 2023. Available from: https://www.bioinformatics.babraham. ac.uk/projects/fastqc/
- 7. Wick RR. 2023. Available from: https://github.com/rrwick/Filtlong
- Wick RR, Judd LM, Gorrie CL, Holt KE. 2017. Unicycler: resolving bacterial genome assemblies from short and long sequencing reads. PLoS Comput Biol 13:e1005595. https://doi.org/10.1371/journal.pcbi.1005595

- Tatusova T, DiCuccio M, Badretdin A, Chetvernin V, Nawrocki EP, Zaslavsky L, Lomsadze A, Pruitt KD, Borodovsky M, Ostell J. 2016. NCBI prokaryotic genome annotation pipeline. Nucleic Acids Res 44:6614– 6624. https://doi.org/10.1093/nar/gkw569
- Meier-Kolthoff JP, Göker M. 2019. TYGS is an automated highthroughput platform for state-of-the-art genome-based taxonomy. Nat Commun 10:2182. https://doi.org/10.1038/s41467-019-10210-3
- Blin K, Shaw S, Augustijn HE, Reitz ZL, Biermann F, Alanjary M, Fetter A, Terlouw BR, Metcalf WW, Helfrich EJN, van Wezel GP, Medema MH, Weber T. 2023. antiSMASH 7.0: new and improved predictions for detection, regulation, chemical structures and visualisation. Nucleic Acids Res 51:W46–W50. https://doi.org/10.1093/nar/gkad344
- Laffont C, Brutesco C, Hajjar C, Cullia G, Fanelli R, Ouerdane L, Cavelier F, Arnoux P. 2019. Simple rules govern the diversity of bacterial nicotianamine-like metallophores. Biochem J 476:2221–2233. https://doi.org/10. 1042/BCJ20190384
- Zavascki AP, Goldani LZ, Li J, Nation RL. 2007. Polymyxin B for the treatment of multidrug-resistant pathogens: a critical review. J Antimicrob Chemother 60:1206–1215. https://doi.org/10.1093/jac/ dkm357