

#### Interaction of the Cysteine-Rich Domain of the Paramecium domesticated transposase PiggyMac with histone H3

Luiza Mamigonian Bessa, Marc Guérineau, Séverine Moriau, Nadine Assrir, Ewen Lescop, François Bontems, Carine van Heijenoort, Éric Guittet, Mireille Bétermier, Julien Bischerour, et al.

#### ▶ To cite this version:

Luiza Mamigonian Bessa, Marc Guérineau, Séverine Moriau, Nadine Assrir, Ewen Lescop, et al.. Interaction of the Cysteine-Rich Domain of the Paramecium domesticated transposase PiggyMac with histone H3. 11th Conference "NMR: a tool for Biology", Jan 2019, Paris, France. 46, pp.2660-2677, 2018, 10.1093/nar/gky044 . hal-04381885

#### HAL Id: hal-04381885 https://cnrs.hal.science/hal-04381885

Submitted on 30 Jan 2024  $\,$ 

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.



# Interaction of the Cysteine-Rich Domain of the *Paramecium* domesticated transposase PiggyMac with histone H3



L. Mamigonian Bessa<sup>\*1</sup>, M. Guérineau<sup>2</sup>, S. Moriau<sup>1</sup>, N. Assrir<sup>1</sup>, E. Lescop<sup>1</sup>, F. Bontems<sup>1</sup>, C. van Heijenoort<sup>1</sup>, E. Guittet<sup>1</sup>, M. Bétermier<sup>2</sup>, J. Bischerour<sup>2</sup>, N. Morellet<sup>1</sup>

<sup>1</sup>Institut de Chimie des Substances Naturelles, CNRS UPR 2301, Université de Paris-Saclay, Gif-Sur-Yvette, France <sup>2</sup>Institute of Integrative Biology of the Cell (I2BC), CEA, CNRS, Université Paris Sud, Université Paris-Saclay, Gif-Sur-Yvette, France



# Background

Transposases catalyze the DNA strand breakage and rejoining elements called DNA transposons. Various organisms have occasionally utilized transposase genes to create new genes with beneficial cellular functions. This is called "domestication" of transposases. Intriguing instances of domestication of piggyBac transposases were found in the ciliate Paramecium tetraurelia.

The PiggyMac (Pgm) domesticated transposase, which carries a potentially active catalytic site (1), is required during the sexual reactions that occur during transposition of mobile genetic cycle of this unicellular eukaryote (Fig. A) for the precise elimination of ≈45,000 short non-coding Internal Eliminated Sequences (IESs) from the polyploid somatic macronucleus (Fig. B). IESs are eliminated through a precise cut-and-close process, during which Pgm cleaves double-stranded DNA at flanking TA dinucleotides located at IES ends.



# A domesticated transposase precisely excises IESs

# Characterization of the Cysteine-rich domain of PiggyMac



## **Objective**

Understand the role of the cysteine-rich domain in PiggyMac's transposase activity.

## **PiggyMac CRD interacts with histone H3 N-terminal tail**





# **PiggyMac CRD does not recognize IES DNA sequences**

IES 1835 (56-mer)

GAATGTTGCGTCGAACGTGTAGAGATCAACTACCATAACACCACAAATAAAATTTTTCATTATATAA-

#### The DNA sequences (dashed lines) were taken from the *P. tetraurelia* genome at or adjacent to excision sites (dark red).



**PiggyMac CRD does not appear to** recognize DNA sequences.

Using two DNA sequences that either contained the entire eliminated sequence (IES 1835) or is upstream (IES 4404) from the excision site (the cleavage site is expected to bind Pgm catalytic domain not the CRD) no changes were observed in the <sup>1</sup>H,<sup>1</sup>H-NOESY spectrum of Pgm CRD.

#### Unlike in PiggyBac (wild transposase from *Trichoplusia ni*) (3), the role of the CRD does not appear to be recognition of the cleavage site by its DNA target sequences





#### Histone H3 interacts with two regions of PiggyMac CRD.

**A**<sup>1</sup>H,<sup>15</sup>N-HSQC spectra (<sup>1</sup>H 800 MHz) of 100 µM <sup>15</sup>N-PiggyMac CRD in the absence (red) and presence of 1 mM histone H3[1-19] (black). **B** Chemical shift perturbations ( $\Delta\Delta\delta$ ) of amide resonances from PiggyMac CRD in spectrum A. C Relative peak intensity of PiggyMac CRD amide resonances, comparing the presence and absence of histone H3. **D** A dissociation constant of 279±70 µM was determined from the chemical shift perturbations. **E** Interaction zones of H3[1-19] shown in blue on PiggyMac CRD structure.

#### ELISA shows interaction between H3 peptides and Pgm CRD.

Biotinylated *P. tetraurelia* histone peptides (see table) were synthesized by Proteogenix. Each well on Pierce Streptavidin Coated High Capacity 96-well Plate was washed 3 times with 200 µL Wash Buffer (25mM Tris pH 7.4; 150mM NaCl; 0.1% BSA; 0.05% Tween 20). The biotinylated peptides (4 µM, in 100 µL Wash Buffer) were applied into wells to saturate their surface and the plate was incubated at room temperature for 1h. The plate was washed three times with 200 µL Wash Buffer and then 100 µL of purified MBP-CRD or MBP (40 nM) were added to each well. The plate was incubated at 4°C overnight before washing (as above) to remove unbound proteins. 100 µL of monoclonal anti-MBP-HRP antibody (NEB #8038S; diluted 1/10000 in Wash Buffer) was added to each well and the plate was incubated for 30 min at room temperature, then washed three times. Peroxidase enzyme activity was measured with tetramethylbenzidine (1 Step<sup>™</sup> Ultra TMB-ELISA, ThermoScientific) using a TECAN Infinite 200 pro

reader at 450 nm.



## An interaction region displays peculiar relaxation behavior



#### NMR relaxation measurements of Pgm CRD.

<sup>15</sup>N R<sub>1</sub>, R<sub>2</sub> and heteronuclear <sup>15</sup>N-<sup>1</sup>H NOE relaxation parameters measured at 950 MHz and 293 K.

The interaction regions identified to interact with H3 are highlighted in blue. The region corresponding to residues 57-71 displays increased R<sub>2</sub>, possibly due to exchange in the µs to ms timescale.

### Conclusions

NMR and ELISA experiments suggest that the CRD may play a role in identifying IES excision sites via the recognition of the histone H3 tail.

Further investigations are needed to understand if there is a preference for the other histones (H1, H2, H4) and/or certain histone modifications.

(1) Dubois E, Mathy N, Régnier V, Bischerour J, Baudry C, Trouslard R and Bétermier M. Multimerization properties of PiggyMac, a domesticated piggyback transposase involved in programmed genome rearrangements. Nucleic Acid Research (2017), 45(6):3204-3216.

- (2) Baudry C, Malinsky S, Restituito M, Kapusta A, Rosa S, Meyer E, Bétermier M. PiggyMac, a domesticated piggyBac transposase involved in programmed genome rearrangements in the ciliate Paramecium tetraurelia. Genes & Development (2009), 23(21):2478-2483
- (3) Morellet N, Li X, Wieninger S A, Taylor J L, Bischerour J, Moriau S, Lescop E, Bardiaux B, Mathy N, Assrir N, Bétermier M, Nilges M, Hickman A B, Dyda F, Craig N L, Guittet E. Sequence-specific DNA binding activity of the cross-brace zinc finger motif of the piggyBac transposase. Nucleic Acids Research (2018) 46(5):2660-2677.